

NYGC Quality Control

All samples undergo rigorous quality assessment using a comprehensive set of quality measures upon completion of each step of sample processing: 1) sample receipt, 2) library preparation, 3) sequencing, 4) data analysis. Quality control measures are scrutinized by the combined efforts of our Project Management, Laboratory, Sequencing Analytics, and Bioinformatics teams. Samples that do not meet our expected quality criteria are flagged and reviewed in consultation with the investigator prior to initiation of the next step of the sample processing pipeline.

RNA Preparation and QC

Extraction

Total RNA is extracted from flash frozen post-mortem tissue. Trizol/Chloroform extraction method is used, followed by Qiagen RNeasy minikit column purification. The column purification step is used to ensure the quality of extracted RNA.

RNA Quantification

Total RNA is quantified using Nanodrop 2000 and Qubit™ 2.0 Fluorometer. Nanodrop reading provides a 260/280 ratio, a generally accepted measure of purity of extracted RNA whereas Qubit reliably measures its concentration.

RNA Integrity

The quality of the RNA preparation and RIN score are verified on an Agilent Bioanalyzer.

Library Preparation and Initial QC

Library Preparation

Libraries are prepared using KAPA Stranded RNA-Seq Kit with RiboErase and unique Illumina-compatible PCR primers with indexes purchased from BioScientific (NEXTflex RNA-seq Barcodes, cat# 512915, 8nt index).

Libraries are prepared using 500ng of total RNA input, 550bp in length and sequenced PE 125, with the yield of ~60 million reads per sample.

Initial Library QC

Libraries are quantified on a Qubit™ 2.0 Fluorometer; the number of amplifiable molecules is estimated by qPCR using the Kapa Library Quantification kit for Illumina Platforms (Kapa Biosystems) and the average size of ~480-550bp verified on an Agilent Bioanalyzer. Libraries with balanced indexes are pooled, subject to qPCR (Kapa Library Quantification kit, Kapa Biosystems) and Qubit quantification before loading onto an Illumina HiSeq 2500 sequencer.

Library QC

Library Quantification

Picogreen is used to measure the total amount of DNA in the prepared library. Quantitative PCR (qPCR) uses specific oligos complimentary to Illumina's TruSeq adapters to measure the amount of adapter-ligated DNA (ligation efficiency) that is compatible with sequencing.

Library Size distribution

Size distribution profiles of the final libraries are assessed using the Fragment analyzer/Bioanalyzer. Libraries that fall outside of the expected size range and/or contain adapter dimer contaminants are flagged.

Sequencing QC

All sequencing runs are reviewed for quality by our Sequencing Analytics team. Sequencing runs that do not pass our quality criteria for each of the metrics below are flagged and reviewed in consultation with Illumina Technical Support.

% Pass Filter (PF) clusters

Library cluster efficiency should fall within the optimal range expected for instrument and flowcell type. PF percentages outside the expected range indicate either incorrect loading concentrations or problems with a particular sequencing run.

% sample de-multiplexed

All PF reads within a single lane of a flowcell are assigned to a specific barcoded library based on the indexed read. The percentage of reads within a lane that are assigned to each sample after de-multiplexing is assessed to confirm expected sample distribution within the sample pool.

of PF reads/sample

The total number of PF reads per sample must meet the expected number of reads for a given sequencing application and analysis type, as discussed upfront with the investigator. Samples that do not meet the expected number of reads are queued for additional sequencing.

% bases >Q30

To ensure the highest quality sequencing data, FASTQ data in which at least 75% (HiSeq X) or 80% (HiSeq 2500) of bases have an Illumina Quality score >30 (a Phred like score indicating an expected 99.9% base call accuracy) are selected and used in downstream analysis.

Quality by cycle

Assessment of the quality score by cycle is used to verify that the accuracy of called bases is maintained across the entire length of the sequencing read.

GC content

GC content is reflective of both sample and library type. GC content can vary between organisms, and can be an indicator of poor sequence quality attributable to biases introduced during library preparation.

K-mer content/adaptor contamination

FASTQC data is examined to identify over-represented sequences in the sequencing data, including k-mers and reads that align to the Illumina adapter sequences, both of which could indicate poor library quality and result in uneven base composition.

Data Analysis QC - RNA

Total reads/sample

The total number of reads per sample must meet the expected total number of reads for a given sequencing application and analysis type, as discussed upfront with the investigator. Samples that do not meet the expected number of reads are queued for additional sequencing.

% rRNA

The proportion of reads in sample that align to ribosomal RNA (rRNA) provides a measure of the success of upfront rRNA (ribo)-depletion during total RNA library preparation, and polyA-enrichment during mRNA library preparation. Higher than expected levels of reads mapping to rRNAs can lower the signal-to-noise ratio, and can have an impact on downstream analysis.

% Duplicates

PCR amplification during library preparation can give rise to the duplication of reads. Libraries that produce a higher than expected number of duplicate reads result from reduced library complexity and reduced representation of the underlying transcript diversity. Higher than expected duplication rates can also indicate reduced levels of sample complexity attributable to lower amounts of starting material used for library preparation.

% Aligned

The proportion of reads successfully mapping to the reference genome.

% Gene assignment

The percentage of mapped reads that are assigned to annotated gene regions are marked as follows: 1) % coding (CDS), 2) % un-translated 3' and 5' regions (UTRs), 3) % intronic (non-coding regions within genes), 4) % intergenic (non-coding outside of genes). Samples that show lower or higher than expected assignment to these genomic regions could indicate quality issues related to technical artifacts introduced during library preparation and/or sequencing, initial quality of the RNA sample, or could reveal information about the origin of the sample (tissue type).

% Strandedness

The percentage of sequenced reads with the correct "strandedness" is verified following preparation of libraries using a stranded library preparation protocol (which preserves DNA strand information).

5'/3' coverage

Read coverage is inspected for uniformity across gene bodies. Ideal experiments show uniform coverage; bias in the representation of reads at either the 5' or 3' ends of gene bodies can indicate poor sample quality and/or technical artifacts introduced during library preparation.

% Mean GC

The mean GC content averaged for each sequenced read is used to flag libraries that show biases introduced through PCR amplification during library preparation and/or sequencing.

Insert size

The mean inner distance between the end of Read 1 and the start of Read2 are calculated to confirm that the library insert sizes are appropriate for the sequenced read length. Smaller insert sizes lead to overlapping reads and/or sequencing into the adapter sequences, limiting the number of usable bases for mapping and downstream analysis.

Xist vs chrY gender check

Gender-specific transcript expression from the X and Y chromosomes is used to determine the gender of the sequenced sample and compare to the gender specified in the sample submission form.

Unsupervised clustering to check for unexpected structure in data (batch effects, sample swap and contamination)

Hierarchical clustering of the RNA sequencing data is performed to validate sample identity and grouping based on experimental design. Unexpected clustering results can reveal potential artifacts in the experiment, such as sample swaps, contamination, sample quality variability and/or technical batch effects, all of which undergo further investigation.